DIAGNOSTIC1525 13th Ave. Seattle, WA 98122CYTOGENETICS800-328-2026fax 206-325-2975INCORPORATEDwww.diagnosticcytogenetics.com

Principal Investigator:	
Department/Company:	
Address:	
Phone:	_Fax:
Submitted By:	
Email:	

SAMPLE INFORMATION

Study Identifier:
Cell Line:
Cell Type:
Passage Date:
Previous Cytogenetic Studies:

TEST REQUESTED

Chromosome Analysis:	Standard 20 cell analysis	Tiered analysis 🛛 Other:	
----------------------	---------------------------	--------------------------	--

□ FISH*____

*Please contact DCI before ordering for more information

BILLING INFORMATION

Please contact DCI for information about billing options.

Billing Contact:	
Phone:	_Fax:
Email:	

SPECIMEN REQUIREMENTS

- 1. Perform a 1:4 split of one confluent T25 flask (approximately 2 million cells per flask).
- 2. The next day, if the cells are adherent and healthy, fill up 2 flasks to capacity with media and send them to us.
- 3. Label the flasks with the cell line, passage date, and your study identifier.
- 4. Please also include an aliquot of media along with the flasks in case it is necessary to feed the culture prior to harvesting for metaphase chromosomes.

SPECIMEN PACKAGING AND SHIPPING

- 1. Before shipping, fill each flask to capacity with warmed media. Tighten the cap and seal it with parafilm.
- 2. Put each flask in a plastic bio-hazard bag with absorbent material.
- 3. Place the specimen bag with the flask in the provided DCI Styrofoam box.
- 4. Include an aliquot of media sealed with parafilm along with the flasks in case it is necessary to feed the culture prior to harvesting for metaphase chromosomes.
- 5. Fill out a DCI Cell Line Services request form and place it in the box.
- 6. Warm ice packs to 37°C and include as extra insulation during shipping.
- 7. Tape the Styrofoam box so that it is securely shut and place in the DCI mailer cardboard box.
- 8. Affix a completed DCI FedEx airbill to the cardboard box for overnight shipment. Please mark all outgoing FedEx samples sent on Friday with a Saturday delivery sticker, as DCI is open for deliveries on Saturday.
- 9. Call FedEx to arrange pick-up at 800-GO-FEDEX (800-463-3339).
- 10. Call DCI with the tracking number.

For Federal Express envelopes, airbills and other supplies, please call DCI at (800) 328-2026.